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References

- Allen WE, Kauvar IV, Chen MZ, Richman EB, Yang SJ, Chan K, Gradinaru V, Deverman BE, Luo L, Deisseroth (2017). Global Representations of Goal-Directed Behavior in Distinct Cell Types of Mouse Neocortex. *Neuron* 94 (4) 891-907, doi.org/10.1016/j.neuron.2017.04.017
- Adams B, Bak H, Tustian AD (2020). Moving from the Bench Towards a Large Scale, Industrial Platform Process for Adeno-Associated Viral Vector Purification. *Biotechnology & Bioengineering*, doi:10.1002/bit.27472
- Chan, K, Jang, M, Yoo, B (2017). Engineered AAVs for efficient noninvasive gene delivery to the central and peripheral nervous systems. *Nat Neurosci* 20, 1172-1179, doi.org/10.1038/nn.4593
- Levy JM, Yeh WH, Pendse N (2020). Cytosine and adenine base editing of the brain, liver, retina, heart and skeletal muscle of mice via adeno-associated viruses. *Nature Biomedical Engineering*, 4(1):97-110, doi:10.1038/s41551-019-0501-5

Q&A

[Q1] How can I inactivate or remove SAN High Quality after use?

[A1] There are various ways of inactivating or removing the protein after use. Due to the high pI of the protein (pI 9.6), it can be removed by the use of cation exchange chromatography. Irreversible inactivation can be achieved by the use of reducing agents (e.g. TCEP, DTT) combined with moderate temperatures. Alternatively, the activity of the enzyme can be inhibited by chelators.

[Q2] Could SAN High Quality be used for purification of viral vectors, e.g. AAV based vectors?

[A2] Yes, this endonuclease is suited for viral vector purifications. The purity of virus vectors is pivotal prior to in vivo transductions of cells or animals. Salt is often an important component to achieve highly purified vectors by minimizing aggregation and increasing target yield. SAN High Quality is highly compatible with the use of high salt conditions.

[Q3] Can I replace my current endonuclease with SAN High Quality in all applications?

[A3] SAN High Quality has a different activity profile and tolerance for certain buffer components compared to other endonucleases. SAN High Quality can replace your current endonuclease in most cases and is a superior nuclease in high salt conditions.

[Q4] What is the advantage of SAN High Quality compared to other nucleases?

[A4] Other nucleases can generally not tolerate salt very well. SAN High Quality has an optimal activity at high salt concentrations, enabling high salt conditions throughout the workflow.

[Q5] What is the lowest reaction temperature I can use for SAN High Quality?

[A5] The optimal reaction temperature of SAN High Quality is 37°C. At lower temperatures, the reaction will be slower. Activity at 4°C is about 5-10% compared to the activity at 37°C. You may compensate by increasing the reaction time or increasing the amount of enzyme.

[Q6] What salt concentrations and pH-values are optimal?

[A6] Optimal conditions of NaCl are between 0.4-0.6 M and optimal pH-values are between 8.5-9.5. Low concentrations of salt can be compensated for by using a higher pH. Low pH can be compensated for by using a higher concentration of salt.

[Q7] Does the product contain BSA?

【A7】 No.

【Q8】 How stable is SAN High Quality?

【A8】 SAN High Quality is a robust enzyme which is stable at room temperature for weeks. The product should be stored at -20°C or 4°C for best storage stability.

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货号	描述	规格	浓度	储存条件
70921-202	SAN HQ (Bioprocessing grade), Triton FREE,25-25k	25 kU	25-30 U/μl	-20 °C
70921-150	SAN HQ (Bioprocessing grade), Triton FREE, 500k	500 kU	>250 U/μl	-20 °C
70921-160	SAN HQ (Bioprocessing grade), Triton FREE, 5M	5 MU	>250 U/μl	-20 °C
70930-001	SAN HQ ELISA Kit	1x96 Well	/	2-8 °C
70950-202	M-SAN HQ (Bioprocessing grade), Triton FREE,25-25k	25 kU	25-30 U/μl	-20 °C
70950-150	M-SAN HQ (Bioprocessing grade), Triton FREE, 500k	500 kU	>250 U/μl	-20 °C
70950-160	M-SAN HQ (Bioprocessing grade), Triton FREE, 5M	5 MU	>250 U/μl	-20 °C
70960-001	M-SAN HQ ELISA Kit	8x12 strips	/	2-8 °C

Cutting-edge enzymes from Norway 来自挪威的尖端酶类

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Shanghai Besto Biotechnology Co., Ltd