

UPULITE® SP C18-PI 着色剂专用柱

微纯Upulite® SP C18-PI着色剂专用柱赋予化合物对称的峰形，柠檬黄与新红、喹啉黄与日落黄、胭脂红等均可基线分离，11种着色剂的保留和分离度满足样品分析定性定量的需求

I 食品中合成着色剂的测定

色谱柱: Upulite® SP C18-PI 5 μ m 4.6*250mm
(货号: SP03-254650)

流速: 0.8mL/min

柱温: 40°C

进样体积: 10 μ L

检测波长: 415nm、520nm、610nm

流动相: A 20mmol/L乙酸铵; B甲醇

梯度:

Time	A%	B%
0	90	10
12	60	40
19	50	50
22.5	45	55
24	5	95
33	5	95
34	90	10
42	90	10

色谱图:

